



E4F1 rabbit pAb

Catalog No	YP-Ab-09242
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	E4F1 E4F
Protein Name	Transcription factor E4F1 (EC 6.3.2.-) (E4F transcription factor 1) (Putative E3 ubiquitin-protein ligase E4F1) (Transcription factor E4F) (p120E4F) (p50E4F)
Immunogen	Synthesized peptide derived from human E4F1
Specificity	This antibody detects endogenous levels of human E4F1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus, nucleoplasm. Cytoplasm. A small fraction is detected in the cytoplasm. Excluded from the nucleolus where it is targeted upon CDKN2A overexpression. Localizes to the mitotic spindle during embryogenesis. .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

