



PDK4 rabbit pAb

| | |
|----------------------------------|---|
| Catalog No | YP-Ab-09350 |
| Isotype | IgG |
| Reactivity | Human;Mouse;Rat |
| Applications | WB |
| Gene Name | PDK4 PDHK4 |
| Protein Name | [Pyruvate dehydrogenase [lipoamide]] kinase isozyme 4, mitochondrial (EC 2.7.11.2) (Pyruvate dehydrogenase kinase isoform 4) |
| Immunogen | Synthesized peptide derived from human PDK4 AA range: 235-285 |
| Specificity | This antibody detects endogenous levels of human PDK4 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Source | Polyclonal, Rabbit,IgG |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Dilution | WB 1:1000-2000 |
| molecular weight | 45kD |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| molecular weight | 45kD |
| Cell Pathway | Mitochondrion matrix. |
| Tissue Specificity | Ubiquitous; highest levels of expression in heart and skeletal muscle. |
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |

Products Images