



# HERC2 rabbit pAb

<b>Catalog No</b>	YP-Ab-09357
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	HERC2
<b>Protein Name</b>	E3 ubiquitin-protein ligase HERC2 (EC 6.3.2.-) (HECT domain and RCC1-like domain-containing protein 2)
<b>Immunogen</b>	Synthesized peptide derived from human HERC2 AA range: 646-696
<b>Specificity</b>	This antibody detects endogenous levels of human HERC2
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	Cytoplasm. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole. Nucleus. Recruited to sites of DNA damage in response to ionizing radiation (IR) via its interaction with RNF8. May loose association with centrosomes during mitosis.
<b>Tissue Specificity</b>	
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images