



# HRSL2 rabbit pAb

<b>Catalog No</b>	YP-Ab-09463
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Rat;Mouse;
<b>Applications</b>	WB
<b>Gene Name</b>	HRASLS2
<b>Protein Name</b>	HRAS-like suppressor 2 (EC 2.3.1.-) (EC 3.1.1.-)
<b>Immunogen</b>	Synthesized peptide derived from human HRSL2 AA range: 1-51
<b>Specificity</b>	This antibody detects endogenous levels of human HRSL2
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	Cytoplasm . Membrane ; Single-pass membrane protein . Exhibits a granular pattern in the cytoplasm with preferential perinuclear localization. .
<b>Tissue Specificity</b>	Expressed in liver, kidney, small intestine testis and colon (PubMed:19615464). Undetectable in testis, placenta, salivary gland and fetal brain (PubMed:18163183).
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

