



MGN2 rabbit pAb

Catalog No	YP-Ab-09661
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MAGOHB MAGOH2
Protein Name	Protein mago nashi homolog 2
Immunogen	Synthesized peptide derived from human MGN2 AA range: 55-105
Specificity	This antibody detects endogenous levels of human MGN2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

