



CPSM rabbit pAb

Catalog No	YP-Ab-09686
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CPS1
Protein Name	Carbamoyl-phosphate synthase [ammonia], mitochondrial (EC 6.3.4.16) (Carbamoyl-phosphate synthetase I) (CPSase I)
Immunogen	Synthesized peptide derived from human CPSM AA range: 816-866
Specificity	This antibody detects endogenous levels of human CPSM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion . Nucleus, nucleolus .
Tissue Specificity	Primarily in the liver and small intestine.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

