



RASN rabbit pAb

Catalog No	YP-Ab-09706
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NRAS HRAS1
Protein Name	GTPase NRas (Transforming protein N-Ras)
Immunogen	Synthesized peptide derived from human RASN AA range: 7-57
Specificity	This antibody detects endogenous levels of human RASN
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane ; Lipid-anchor ; Cytoplasmic side . Golgi apparatus membrane ; Lipid-anchor . Shuttles between the plasma membrane and the Golgi apparatus. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

