



# TRAT1 rabbit pAb

<b>Catalog No</b>	YP-Ab-09708
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	TRAT1 TCRIM HSPC062
<b>Protein Name</b>	T-cell receptor-associated transmembrane adapter 1 (T-cell receptor-interacting molecule) (TRIM) (pp29/30)
<b>Immunogen</b>	Synthesized peptide derived from human TRAT1 AA range: 56-106
<b>Specificity</b>	This antibody detects endogenous levels of human TRAT1
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	Cell membrane ; Single-pass type III membrane protein .
<b>Tissue Specificity</b>	Strongly expressed in thymus, and to a lesser extent in spleen, lymph node and peripheral blood lymphocytes. Present in T-cells and NK cells, but not B-cells (at protein level).
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

