



DJB14 rabbit pAb

| | |
|----------------------------------|--|
| Catalog No | YP-Ab-09728 |
| Isotype | IgG |
| Reactivity | Human;Mouse |
| Applications | WB |
| Gene Name | DNAJB14 UNQ9427/PRO34683 |
| Protein Name | DnaJ homolog subfamily B member 14 |
| Immunogen | Synthesized peptide derived from human DJB14 AA range: 159-209 |
| Specificity | This antibody detects endogenous levels of human DJB14 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Source | Polyclonal, Rabbit,IgG |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Dilution | WB 1:1000-2000 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | |
| Observed Band | |
| Cell Pathway | Endoplasmic reticulum membrane ; Single-pass membrane protein . Nucleus membrane ; Single-pass membrane protein . Localizes to the endoplasmic reticulum membrane (PubMed:23018488, PubMed:24732912, PubMed:27916661). When overexpressed, forms membranous structures in the nucleus (PubMed:24732912). . |
| Tissue Specificity | |
| Function | |
| Background | |
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



Products Images