



CF089 rabbit pAb

Catalog No	YP-Ab-09736
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C6orf89 BRAP UNQ177/PRO203
Protein Name	Bombesin receptor-activated protein C6orf89 (Amfion)
Immunogen	Synthesized peptide derived from human CF089 AA range: 193-243
Specificity	This antibody detects endogenous levels of human CF089
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	[Isoform 1]: Golgi apparatus membrane ; Single-pass type II membrane protein . Midbody .; [Isoform 2]: Golgi apparatus membrane ; Single-pass type II membrane protein . Midbody .; [Isoform 3]: Cytoplasm . Nucleus, nucleolus . Retained in nucleolar organiser regions (NORs) in mitotic cells (PubMed:23460338). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images