



CLU rabbit pAb

Catalog No	YP-Ab-09738
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KIAA0664
Protein Name	Clustered mitochondria protein homolog
Immunogen	Synthesized peptide derived from human CLU AA range: 435-485
Specificity	This antibody detects endogenous levels of human CLU
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
molecular weight	144kD
Cell Pathway	Cytoplasm . Cytoplasmic granule . A fraction colocalizes with tyrosinated tubulin and can be detected close to mitochondria. .
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images