



LNx1 rabbit pAb

Catalog No	YP-Ab-09753
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LNx1 LNx PDZRN2 UNQ574/PRO1136
Protein Name	E3 ubiquitin-protein ligase LNx (EC 6.3.2.-) (Ligand of Numb-protein X 1) (Numb-binding protein 1) (PDZ domain-containing RING finger protein 2)
Immunogen	Synthesized peptide derived from human LNx1 AA range: 527-577
Specificity	This antibody detects endogenous levels of human LNx1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed in heart, placenta, kidney, pancreas and brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

