



# NCKPL rabbit pAb

<b>Catalog No</b>	YP-Ab-09762
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Rat;Mouse;
<b>Applications</b>	WB
<b>Gene Name</b>	NCKAP1L HEM1
<b>Protein Name</b>	Nck-associated protein 1-like (Hematopoietic protein 1) (Membrane-associated protein HEM-1)
<b>Immunogen</b>	Synthesized peptide derived from human NCKPL AA range: 582-632
<b>Specificity</b>	This antibody detects endogenous levels of human NCKPL
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	Cell membrane ; Single-pass membrane protein ; Cytoplasmic side . Cytoplasm . Localizes to the leading edge of polarized neutrophils. .
<b>Tissue Specificity</b>	Expressed only in cells of hematopoietic origin (PubMed:7643388, PubMed:1932118). Expressed in neutrophils (at protein level) (PubMed:16417406). Expressed in T-cells (at protein level) (PubMed:32647003).
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images