



MP17L rabbit pAb

Catalog No	YP-Ab-09773
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MPV17L
Protein Name	Mpv17-like protein (M-LP homolog) (M-LPH)
Immunogen	Synthesized peptide derived from human MP17L AA range: 39-89
Specificity	This antibody detects endogenous levels of human MP17L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	[Isoform 1]: Peroxisome membrane ; Multi-pass membrane protein.
Tissue Specificity	Isoform 1 is detected in the kidney (at protein level). Isoform 1 and isoform 2 are expressed in the kidney, heart, liver, lung, pancreas and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

