



kAZD1 rabbit pAb

Catalog No	YP-Ab-09807
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KAZALD1 FKSG28 FKSG40 UNQ2945/PRO21184
Protein Name	Kazal-type serine protease inhibitor domain-containing protein 1
Immunogen	Synthesized peptide derived from human KAZD1 AA range: 4-54
Specificity	This antibody detects endogenous levels of human KAZD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

