



TXLNB rabbit pAb

Catalog No	YP-Ab-09830
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TXLNB C6orf198 MDP77
Protein Name	Beta-taxilin (Muscle-derived protein 77) (hMDP77)
Immunogen	Synthesized peptide derived from human TXLNB AA range: 425-475
Specificity	This antibody detects endogenous levels of human TXLNB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Expressed in skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

