



LIPJ rabbit pAb

Catalog No	YP-Ab-09853
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LIPJ LIPL1
Protein Name	Lipase member J (EC 3.1.1.-) (Lipase-like abhydrolase domain-containing protein 1)
Immunogen	Synthesized peptide derived from human LIPJ AA range: 228-278
Specificity	This antibody detects endogenous levels of human LIPJ
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

