



CV032 rabbit pAb

Catalog No	YP-Ab-09905
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	C22orf32
Protein Name	UPF0466 protein C22orf32, mitochondrial
Immunogen	Synthesized peptide derived from human CV032 AA range: 9-59
Specificity	This antibody detects endogenous levels of human CV032
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion inner membrane ; Single-pass membrane protein . MAIP1 is required to assist sorting of EMRE/SMDT1 into mitochondrion by protecting EMRE/SMDT1 against protein degradation by YME1L1, thereby ensuring SMDT1/EMRE maturation by the mitochondrial processing peptidase (PMPCA and PMPCB) (PubMed:27642048). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images