



HXD4 rabbit pAb

Catalog No	YP-Ab-09923
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	HOXD4 HOX4B
Protein Name	Homeobox protein Hox-D4 (Homeobox protein HHO.C13) (Homeobox protein Hox-4B) (Homeobox protein Hox-5.1)
Immunogen	Synthesized peptide derived from human HXD4 AA range: 99-149
Specificity	This antibody detects endogenous levels of human HXD4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

