



UVSSA rabbit pAb

Catalog No	YP-Ab-09949
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UVSSA KIAA1530
Protein Name	UV-stimulated scaffold protein A
Immunogen	Synthesized peptide derived from human UVSSA AA range: 516-566
Specificity	This antibody detects endogenous levels of human UVSSA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Chromosome . Accumulates at UV DNA damage sites. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

