



TROAP rabbit pAb

Catalog No	YP-Ab-09960
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TROAP
Protein Name	Tastin (Trophinin-assisting protein) (Trophinin-associated protein)
Immunogen	Synthesized peptide derived from human TROAP AA range: 532-582
Specificity	This antibody detects endogenous levels of human TROAP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm.
Tissue Specificity	Strong expression at implantation sites. Was exclusively localized to the apical side of the syncytiotrophoblast. Also found in macrophages.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

