



# GOG6A rabbit pAb

<b>Catalog No</b>	YP-Ab-09969
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Rat;Mouse;
<b>Applications</b>	WB
<b>Gene Name</b>	GOLGA6A GLP GOLGA6
<b>Protein Name</b>	Golgin subfamily A member 6A (Golgin linked to PML) (Golgin-like protein)
<b>Immunogen</b>	Synthesized peptide derived from human GOG6A AA range: 123-173
<b>Specificity</b>	This antibody detects endogenous levels of human GOG6A
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	
<b>Tissue Specificity</b>	Highly expressed in seminiferous tubes in testis. Highly expressed in spermatids, barely detectable in late pachytene spermatocytes, and not detectable in spermatogonia. Detected at intermediate levels in pancreas and lymph nodes, and at much lower levels in spleen, peripheral blood leukocytes, skeletal muscle, liver, lung, placenta, brain and heart.
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images