



HEM0 rabbit pAb

Catalog No	YP-Ab-09986
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ALAS2 ALASE ASB
Protein Name	5-aminolevulinate synthase, erythroid-specific, mitochondrial (ALAS-E) (EC 2.3.1.37) (5-aminolevulinic acid synthase 2) (Delta-ALA synthase 2) (Delta-aminolevulinate synthase 2)
Immunogen	Synthesized peptide derived from human HEM0 AA range: 434-484
Specificity	This antibody detects endogenous levels of human HEM0
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion inner membrane ; Peripheral membrane protein . Localizes to the matrix side of the mitochondrion inner membrane. . ; [Isoform 4]: Mitochondrion inner membrane ; Peripheral membrane protein . Localizes to the matrix side of the mitochondrion inner membrane. .
Tissue Specificity	[Isoform 1]: Erythroid-specific. ; [Isoform 2]: Erythroid-specific. ; [Isoform 3]: Erythroid-specific. ; [Isoform 4]: Erythroid-specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images