



RDH8 rabbit pAb

Catalog No	YP-Ab-09989
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	RDH8 PRRDH
Protein Name	Retinol dehydrogenase 8 (EC 1.1.1.300) (Photoreceptor outer segment all-trans retinol dehydrogenase)
Immunogen	Synthesized peptide derived from human RDH8 AA range: 14-64
Specificity	This antibody detects endogenous levels of human RDH8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Detected in photoreceptor outer segments in the retina (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

