



VAC14 rabbit pAb

Catalog No	YP-Ab-10037
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VAC14 TAX1BP2 TRX
Protein Name	Protein VAC14 homolog (Tax1-binding protein 2)
Immunogen	Synthesized peptide derived from human VAC14 AA range: 569-619
Specificity	This antibody detects endogenous levels of human VAC14
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endosome membrane . Microsome membrane . Mainly associated with membranes of the late endocytic pathway.
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

