



# KLH31 rabbit pAb

<b>Catalog No</b>	YP-Ab-10066
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	KLHL31 BKLHD6 KBTBD1 KLHL
<b>Protein Name</b>	Kelch-like protein 31 (BTB and kelch domain-containing protein 6) (Kelch repeat and BTB domain-containing protein 1) (Kelch-like protein KLHL)
<b>Immunogen</b>	Synthesized peptide derived from human KLH31 AA range: 342-392
<b>Specificity</b>	This antibody detects endogenous levels of human KLH31
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	
<b>Tissue Specificity</b>	Strongly expressed in skeletal muscle and weakly in heart. According to PubMed:15302408, not expressed in other tissues. According to PubMed:18719355, abundantly expressed in both embryonic skeletal and heart tissues.
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images