



PNMA2 rabbit pAb

Catalog No	YP-Ab-10134
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PNMA2 KIAA0883 MA2
Protein Name	Paraneoplastic antigen Ma2 (40 kDa neuronal protein) (Onconeural antigen Ma2) (Paraneoplastic neuronal antigen MM2)
Immunogen	Synthesized peptide derived from human PNMA2 AA range: 153-203
Specificity	This antibody detects endogenous levels of human PNMA2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Brain-specific. In some cancer patients, specifically expressed by testicular tumor cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

