



AP4M1 rabbit pAb

Catalog No	YP-Ab-10166
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AP4M1 MUARP2
Protein Name	AP-4 complex subunit mu-1 (AP-4 adapter complex mu subunit) (Adapter-related protein complex 4 mu-1 subunit) (Mu subunit of AP-4) (Mu-adaptin-related protein 2) (mu-ARP2) (Mu4-adaptin) (mu4)
Immunogen	Synthesized peptide derived from human AP4M1 AA range: 312-362
Specificity	This antibody detects endogenous levels of human AP4M1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Golgi apparatus, trans-Golgi network membrane ; Peripheral membrane protein . Early endosome . Found in soma and dendritic shafts of neuronal cells. .
Tissue Specificity	Ubiquitous. Highly expressed in testis and lowly expressed in brain and lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

