



HIG1B rabbit pAb

Catalog No	YP-Ab-10174
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	HIGD1B
Protein Name	HIG1 domain family member 1B (Protein CLST 11240)
Immunogen	Synthesized peptide derived from human HIG1B AA range: 19-69
Specificity	This antibody detects endogenous levels of human HIG1B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
molecular weight	11kD
Observed Band	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

