



YPEL3 rabbit pAb

| | |
|----------------------------------|---|
| Catalog No | YP-Ab-10175 |
| Isotype | IgG |
| Reactivity | Human;Mouse |
| Applications | IHC |
| Gene Name | YPEL3 FKSG5 |
| Protein Name | Protein yippee-like 3 |
| Immunogen | Synthesized peptide derived from human YPEL3. AA range:100-119 |
| Specificity | This antibody detects endogenous levels of human YPEL3 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Source | Polyclonal, Rabbit,IgG |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Dilution | IHC-p 1: 50-200 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | |
| Observed Band | |
| Cell Pathway | Nucleus, nucleolus . |
| Tissue Specificity | Widely expressed. |
| Function | |
| Background | |
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



Products Images

