



CPEB3 rabbit pAb

Catalog No	YP-Ab-10188
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CPEB3 KIAA0940
Protein Name	Cytoplasmic polyadenylation element-binding protein 3 (CPE-BP3) (CPE-binding protein 3) (hCPEB-3)
Immunogen	Synthesized peptide derived from human CPEB3 AA range: 443-493
Specificity	This antibody detects endogenous levels of human CPEB3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Nucleus . Cell junction, synapse . Cell projection, dendrite . Cell junction, synapse, postsynaptic density . Predominantly cytoplasmic in unstimulated neurons but translocates to the nucleus following neuronal stimulation (PubMed:20639532, PubMed:22730302). Nuclear import is mediated by importin IPO5 (By similarity). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images