



ERO1B rabbit pAb

Catalog No	YP-Ab-10191
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ERO1LB
Protein Name	ERO1-like protein beta (ERO1-L-beta) (EC 1.8.4.-) (Endoplasmic oxidoreductin-1-like protein B) (Oxidoreductin-1-L-beta)
Immunogen	Synthesized peptide derived from human ERO1B AA range: 23-73
Specificity	This antibody detects endogenous levels of human ERO1B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endoplasmic reticulum membrane ; Peripheral membrane protein ; Luminal side . The association with ERP44 may be essential for its retention in the endoplasmic reticulum.
Tissue Specificity	Highly expressed in the digestive tract, including the duodenum and lower digestive tract. In the stomach, highly expressed in enzyme-producing chief cells (at protein level). In the pancreas, expressed in islets of Langerhans and, at lower levels, in enzyme-secreting cells (at protein level). Detected at low level in many other tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images