



GID8 rabbit pAb

Catalog No	YP-Ab-10203
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	GID8 C20orf11
Protein Name	Glucose-induced degradation protein 8 homolog (Two hybrid-associated protein 1 with RanBPM) (Twa1)
Immunogen	Synthesized peptide derived from human GID8 AA range: 115-165
Specificity	This antibody detects endogenous levels of human GID8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Nucleus . Localizes in the cytoplasm in the absence of Wnt stimulation and in the nucleus in the presence of Wnt stimulation. .
Tissue Specificity	Up-regulated in colorectal cancer tissues (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

