



AP4S1 rabbit pAb

Catalog No	YP-Ab-10209
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	AP4S1
Protein Name	AP-4 complex subunit sigma-1 (AP-4 adapter complex subunit sigma-1) (Adapter-related protein complex 4 subunit sigma-1) (Sigma-1 subunit of AP-4) (Sigma-4-adaptin) (Sigma4-adaptin)
Immunogen	Synthesized peptide derived from human AP4S1 AA range: 69-119
Specificity	This antibody detects endogenous levels of human AP4S1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Golgi apparatus, trans-Golgi network membrane ; Peripheral membrane protein .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

