



NAGK (phospho-Tyr205) rabbit pAb

Catalog No	YP-Ab-10432
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	NAGK
Protein Name	NAGK (phospho-Tyr205)
Immunogen	Synthesized peptide derived from human NAGK (phospho-Tyr205)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat NAGK (phospho-Tyr205)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	N-acetyl-D-glucosamine kinase (N-acetylglucosamine kinase;EC 2.7.1.59;GlcNAc kinase)
Observed Band	
Cell Pathway	
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images