




BCLAF1 (Phospho-Ser531) rabbit pAb

Catalog No	YP-Ab-10552
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	BCLAF1 BTF KIAA0164
Protein Name	BCLAF1 (Phospho-Ser531)
Immunogen	Synthesized peptide derived from human BCLAF1 (Phospho-Ser531)
Specificity	This antibody detects endogenous levels of BCLAF1 (Phospho-Ser531) at Human, Mouse, Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.208% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Bcl-2-associated transcription factor 1 (Btf)
Observed Band	
Cell Pathway	Cytoplasm. Nucleus. Nucleus speckle . Nucleus, nucleoplasm .
Tissue Specificity	Ubiquitous.
Function	function: Death-promoting transcriptional repressor., sequence caution: Contaminating sequence. Potential poly-A sequence., subunit: Interacts with Bcl-2 related proteins, EMD, with the adenovirus E1B 19 kDa protein and with DNA., tissue specificity: Ubiquitous.,
Background	This gene encodes a transcriptional repressor that interacts with several members of the BCL2 family of proteins. Overexpression of this protein induces apoptosis, which can be suppressed by co-expression of BCL2 proteins. The protein localizes to dot-like structures throughout the nucleus, and redistributes to a zone near the nuclear envelope in cells undergoing apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various,using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000