



PLEKHG5 (Phospho-ser881) rabbit pAb

Catalog No	YP-Ab-10566
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	PLEKHG5 KIAA0720
Protein Name	PLEKHG5 (Phospho-ser881)
Immunogen	Synthesized peptide derived from human PLEKHG5 (Phospho-ser881)
Specificity	This antibody detects endogenous levels of PLEKHG5 (Phospho-ser881) at Human, Mouse, Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.222% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Pleckstrin homology domain-containing family G member 5 (PH domain-containing family G member 5) (Guanine nucleotide exchange factor 720) (GEF720)
Observed Band	
Cell Pathway	Cytoplasm . Cytoplasm, perinuclear region . Cell membrane . Cell junction . Cell projection, lamellipodium . Predominantly cytoplasmic, however when endothelial cells are stimulated with lysophosphatidic acid, PLEKHG5 is found in perinuclear regions and at the cell membrane. Localizes at cell-cell junctions in quiescent endothelial cells, and relocalizes to cytoplasmic vesicle and the leading edge of lamellipodia in migrating endothelial cells. .
Tissue Specificity	Predominantly expressed in the peripheral nervous system and brain. Highest expression is observed in heart, lung, kidney, testis and moderate expression is present in spleen, pancreas, skeletal muscle, ovary and liver. Weakly expressed in glioblastoma (GBM) cell lines.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images