



Probable ganciclovir kinase U69 Polyclona Antibody

Catalog No	YP-Ab-10883
Isotype	IgG
Reactivity	Human virus
Applications	ELISA
Gene Name	U69
Protein Name	Probable ganciclovir kinase U69
Immunogen	Synthesized peptide derived from human Probable ganciclovir kinase U69
Specificity	This antibody detects endogenous levels of human Probable ganciclovir kinase U69
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	ELISA(peptide) 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	P52446 Phosphotransferase(U69) Human herpesvirus 6B function:Phosphorylates the antiviral nucleoside analog ganciclovir.,similarity:Belongs to the protein kinase superfamily. Tyr protein kinase family. HCMV ganciclovir subfamily.,
Background	P52446 Phosphotransferase(U69) Human herpesvirus 6B function:Phosphorylates the antiviral nucleoside analog ganciclovir.,similarity:Belongs to the protein kinase superfamily. Tyr protein kinase family. HCMV ganciclovir subfamily.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images