




TMA16 rabbit pAb

Catalog No	YP-Ab-10938
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMA16 C4orf43
Protein Name	TMA16
Immunogen	Synthesized peptide derived from human TMA16 AA range: 147-197
Specificity	This antibody detects endogenous levels of TMA16 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	nucleus,nucleolus,ribosome,
Tissue Specificity	
Function	similarity:Belongs to the UPF0534 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Jarkat cells,
primary antibody was diluted at 1:1000, 4°over night