




TM201 rabbit pAb

Catalog No	YP-Ab-11262
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMEM201 NET5 SAMP1
Protein Name	TM201
Immunogen	Synthesized peptide derived from human TM201 AA range: 290-340
Specificity	This antibody detects endogenous levels of TM201 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	[Isoform 2]: Nucleus inner membrane ; Multi-pass membrane protein . Cytoplasm, cytoskeleton, spindle pole . The C-terminal of isoform 2 is located on the nucleoplasmic side. During interphase, isoform 2 is distributed in the inner nuclear membrane in distinct micro-domains and during mitosis, it is found in the ER but it also localizes to the polar regions of the mitotic spindle. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night