




LITD1 rabbit pAb


Catalog No	YP-Ab-11307
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;IHC
Gene Name	L1TD1 ECAT11
Protein Name	LITD1
Immunogen	Synthesized peptide derived from human LITD1 AA range: 797-847
Specificity	This antibody detects endogenous levels of LITD1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	intracellular ribonucleoprotein complex,
Tissue Specificity	
Function	sequence caution:Translated as Leu.,sequence caution:Translation N-terminally extended.,similarity:Belongs to the transposase 22 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human meningioma. 1, Antibody was diluted at 1:200(4 ° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).