




FA81B rabbit pAb

Catalog No	YP-Ab-11314
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FAM81B
Protein Name	FA81B
Immunogen	Synthesized peptide derived from human FA81B AA range: 33-83
Specificity	This antibody detects endogenous levels of FA81B at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	nucleus,
Tissue Specificity	
Function	caution:It is uncertain whether Met-1 or Met-39 is the initiator. Met-1 does not seem to exist in orthologs.,similarity:Belongs to the FAM81 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4° over night