




C2D2B rabbit pAb

Catalog No	YP-Ab-11341
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CC2D2B C10orf130
Protein Name	C2D2B
Immunogen	Synthesized peptide derived from human C2D2B AA range: 81-131
Specificity	This antibody detects endogenous levels of C2D2B at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4°over night