




DET1 rabbit pAb

Catalog No	YP-Ab-11619
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	DET1
Protein Name	DET1
Immunogen	Synthesized peptide derived from human DET1 AA range: 48-98
Specificity	This antibody detects endogenous levels of DET1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in the ovary, some lymphoid organs and resting leukocytes.
Function	function:Component of the E3 ubiquitin ligase DCX DET1-COP1 complex, which is required for ubiquitination and subsequent degradation of target proteins. The complex is involved in JUN ubiquitination and degradation.,similarity:Belongs to the DET1 family.,subunit:Component of the DCX DET1-COP1 ubiquitin ligase complex at least composed of RBX1, DET1, DDB1, CUL4A and COP1.,tissue specificity:Highly expressed in the ovary, some lymphoid organs and resting leukocytes.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night