




AGAP4 rabbit pAb

Catalog No	YP-Ab-11676
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	AGAP4 CTGLF1 MRIP2
Protein Name	AGAP4
Immunogen	Synthesized peptide derived from human AGAP4 AA range: 296-346
Specificity	This antibody detects endogenous levels of AGAP4 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	function:Putative GTPase-activating protein .,miscellaneous:Encoded by one of the numerous copies of centaurin gamma-like genes clustered in the q11 region of chromosome 10.,similarity:Belongs to the centaurin gamma-like family.,similarity:Contains 1 Arf-GAP domain.,similarity:Contains 1 PH domain.,similarity:Contains 2 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4° over night