




WDR49 rabbit pAb

Catalog No	YP-Ab-11779
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	WDR49
Protein Name	WDR49
Immunogen	Synthesized peptide derived from human WDR49 AA range: 345-395
Specificity	This antibody detects endogenous levels of WDR49 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 8 WD repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SW480 cells,
primary antibody was diluted at 1:1000, 4° over night