




AGAP8 rabbit pAb

Catalog No	YP-Ab-11825
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	AGAP8 CTGLF5
Protein Name	AGAP8
Immunogen	Synthesized peptide derived from human AGAP8 AA range: 339-389
Specificity	This antibody detects endogenous levels of AGAP8 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night