




# DJC21 rabbit pAb

|                                  |                                                                                                                                                                                       |
|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                | YP-Ab-12045                                                                                                                                                                           |
| <b>Isotype</b>                   | IgG                                                                                                                                                                                   |
| <b>Reactivity</b>                | Human;Rat;Mouse;                                                                                                                                                                      |
| <b>Applications</b>              | WB                                                                                                                                                                                    |
| <b>Gene Name</b>                 | DNAJC21 DNAJA5                                                                                                                                                                        |
| <b>Protein Name</b>              | DJC21                                                                                                                                                                                 |
| <b>Immunogen</b>                 | Synthesized peptide derived from human DJC21 AA range: 255-305                                                                                                                        |
| <b>Specificity</b>               | This antibody detects endogenous levels of DJC21 at Human                                                                                                                             |
| <b>Formulation</b>               | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                               |
| <b>Source</b>                    | Polyclonal, Rabbit,IgG                                                                                                                                                                |
| <b>Purification</b>              | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                                             |
| <b>Dilution</b>                  | WB 1: 500-2000                                                                                                                                                                        |
| <b>Concentration</b>             | 1 mg/ml                                                                                                                                                                               |
| <b>Purity</b>                    | ≥90%                                                                                                                                                                                  |
| <b>Storage Stability</b>         | -20°C/1 year                                                                                                                                                                          |
| <b>Synonyms</b>                  |                                                                                                                                                                                       |
| <b>Observed Band</b>             |                                                                                                                                                                                       |
| <b>Cell Pathway</b>              | Cytoplasm . Nucleus . Nucleus, nucleolus . Within the nucleus, localizes primarily to the nucleolus. .                                                                                |
| <b>Tissue Specificity</b>        | Expressed in brain, placenta, kidney and pancreas.                                                                                                                                    |
| <b>Function</b>                  | function:May act as a co-chaperone for HSP70.,similarity:Contains 1 J domain.,similarity:Contains 2 C2H2-type zinc fingers.,tissue specificity:Brain, placenta, kidney and pancreas., |
| <b>Background</b>                |                                                                                                                                                                                       |
| <b>matters needing attention</b> | Avoid repeated freezing and thawing!                                                                                                                                                  |
| <b>Usage suggestions</b>         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                                                     |



## Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night