




# PPR32 rabbit pAb


|                                  |                                                                                                                                                                                                                                            |
|----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                | YP-Ab-12132                                                                                                                                                                                                                                |
| <b>Isotype</b>                   | IgG                                                                                                                                                                                                                                        |
| <b>Reactivity</b>                | Human; Mouse; Rat                                                                                                                                                                                                                          |
| <b>Applications</b>              | WB; IHC                                                                                                                                                                                                                                    |
| <b>Gene Name</b>                 | PPP1R32 C11orf66                                                                                                                                                                                                                           |
| <b>Protein Name</b>              | PPR32                                                                                                                                                                                                                                      |
| <b>Immunogen</b>                 | Synthesized peptide derived from human PPR32 AA range: 72-122                                                                                                                                                                              |
| <b>Specificity</b>               | This antibody detects endogenous levels of PPR32 at Human/Mouse/Rat                                                                                                                                                                        |
| <b>Formulation</b>               | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                    |
| <b>Source</b>                    | Polyclonal, Rabbit, IgG                                                                                                                                                                                                                    |
| <b>Purification</b>              | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                                                                                                  |
| <b>Dilution</b>                  | WB 1:500-2000; IHC-p 1:50-300                                                                                                                                                                                                              |
| <b>Concentration</b>             | 1 mg/ml                                                                                                                                                                                                                                    |
| <b>Purity</b>                    | ≥90%                                                                                                                                                                                                                                       |
| <b>Storage Stability</b>         | -20°C/1 year                                                                                                                                                                                                                               |
| <b>Synonyms</b>                  |                                                                                                                                                                                                                                            |
| <b>Observed Band</b>             |                                                                                                                                                                                                                                            |
| <b>Cell Pathway</b>              | Cell projection, cilium . Cytoplasm . Localized to the cilia of polarized ependymal cells during development and at the adult stage. Preferentially locates between the peripheral microtubules of the axoneme and the ciliary membrane. . |
| <b>Tissue Specificity</b>        |                                                                                                                                                                                                                                            |
| <b>Function</b>                  |                                                                                                                                                                                                                                            |
| <b>Background</b>                |                                                                                                                                                                                                                                            |
| <b>matters needing attention</b> | Avoid repeated freezing and thawing!                                                                                                                                                                                                       |
| <b>Usage suggestions</b>         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                                                                                                          |



## Products Images



Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight). 2, Tris-EDTA, pH 9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 45 min).