



YPEL3 mouse mAb

Catalog No	YP-mAb-10175
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC
Gene Name	YPEL3 FKSG5
Protein Name	Protein yippee-like 3
Immunogen	Synthesized peptide derived from human YPEL3. AA range:100-119
Specificity	This antibody detects endogenous levels of human YPEL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1: 50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

